

RELIABILITY TEST REPORT

Report No: 20201006001

Page: 1 of 4

Date: October 6, 2020

RAYTAC CORPORATION
5F., NO.3, JIANKANG ROAD, ZHONGHE DISTRICT,
NEW TAIPEI CITY, 23586, TAIWAN

Product Description: nRF52833 Bluetooth Module
Style / Item No.: MDBT50-512K/P512K / No.1, No.2
Manufacturer/ Vendor: Raytac Corporation
Country of Origin: Taiwan
Testing Required: High temperature maintain test
Testing Method: JASO D001-94, Clause 5.14
Quantity: Total 2 sets
Testing Period: Sep. 14, 2020 to Sep. 17, 2020
Testing Results: -Please See Attached Sheets

**Signed for and on behalf of
Raytac Corporation.**

Stanley Huang
Deputy Manager

RELIABILITY TEST REPORT

Report No: 20201006001

Page: 2 of 4

Date: October 6, 2020

1. High temperature maintain test:

Test Equipment:

Name	Brands	Model	Serial No.
Testing Chamber	Haich	HC-TOPH-22-CH	QHC20200707-122-3
DC Power Supply	GOODWILL	GPC-1850D	E873364

Lab Environmental Conditions:

Ambient Temperature: $(25 \pm 3) ^\circ\text{C}$

Ambient humidity: $(55 \pm 20) \% \text{ RH}$

Test Method/ Specification:

Test method: Refer to JASO D001-94, Clause 5.14

Sample condition: Operating

Temperature: 105°C

Duration: 72 Hours

- Sample Condition: Operating
- Examine the appearance of specimen(s) by visual check and perform functional check during this test.
- Functional check: Connect the specimen(s) with notebook and examine whether the signal transmission could work normally or not.



RELIABILITY TEST REPORT

Report No: 20201006001

Page: 3 of 4

Date: October 6, 2020

Specimen:

Style / Item No.: MDBT50-512K/P512K / No.1, No.2

Quantity: Total 2 sets

Testing time: 72 hours

Data packet size: 244 bytes

Data rate: 244 bytes * 4 packets/sec = 976 bytes/sec

Number of data sent: Total 252,979,200 bytes

Test Result:

Check Item Style/Items No.	Appearance check (Visual check)	Functional Check
MDBT50-512K / No.1	No visible damage	No data loss
MDBT50-P512K / No.2	No visible damage	No data loss


RELIABILITY TEST REPORT

Report No: 20201006001

Page: 4 of 4

Date: October 6,2020


Test Photos:




1. Appearance of specimen




2. High temperature maintain test



3. High temperature maintain test



4. Functional check
MDBT50-512K / No.1



5. Functional check
MDBT50-P512K / No.2

The End of Test Report